



Scuba Spec, Inc.
220 Richland Ave., W.
Aiken, SC 29801
1-800-457-5155
Fax: 803-643-8056

SCUBA SPEC
199EB

The **Scuba Spec 199EB** is an underwater spectacle kit designed specifically for the visually deficient diver. Prescription lenses are edged and inserted into the frame and snapped in and out of the mask at the diver's discretion. The same prescription (Rx) power as in the diver's habitually worn glasses (including bifocals, if desired) is used. There is no vertex variation.

INSTALLATION INSTRUCTIONS

The anchor bracket is bonded to the face plate using a silicon adhesive cement (e.g., Dow Corning aquarium sealer). To install the anchor bracket, the following procedure is used:

- Thoroughly clean and dry the surface area of the faceplate where the bracket will be installed. It is recommended that you use water or a mild soapy solution. Do not use window cleaner or an alcohol-based solution.
- Clean and roughen the anchor bracket's bonding surface (e.g., with fine sandpaper). Make sure that the four holes in the anchor bracket are clear (unclogged).
- Preposition the anchor bracket on the inner surface of the faceplate so that you know where to bond it. Proper position for the anchor bracket is just above the mask's 'nose pocket' (Figure 1).
- Squeeze a generous amount of silicon adhesive cement onto the back of the anchor bracket.
- Press the anchor bracket onto the faceplate, firmly applying pressure for approximately one minute. Make sure that the anchor bracket is centered and level and that all air bubbles have been expressed from between the bonding surfaces. A small amount of adhesive should extrude through the holes in the anchor bracket. This will help to create a superior lock on the faceplate.
- Allow at least 24 hours for curing of the adhesive bond.
- Use a single-edged razor blade or other suitable tool to remove excess adhesive from around the anchor bracket. Also, remove any dried adhesive protruding from the anchor bracket holes so that the frame front will snap securely into place.



Figure 1. Proper Bracket Location

To unsnap the frame from the anchor bracket, grasp each eye rim between the thumb and middle finger and gently rock the frame off of the bracket.

To prevent fogging of the face plate and the Rx lenses, a good anti-fog preparation should be applied before diving or snorkeling. Inquire about same from your eye care professional or dive shop. An alternative surfactant such as Lemon Joy detergent can also be used. Apply the detergent to the glass surfaces and briskly 'polish' it on until the glass surfaces are 'squeaky-clean'.

NOTE TO OPTICIAN:

- The "R" and "L" eye designations imprinted on the back (flat side) of the frame is a mold error. They are reversed. Please edge accordingly.
- The nylon frame front should be hydrated every 2-3 months by soaking in water for 24 hours. This will help maintain the frame's flexibility.

FRAME FRONT AND BRACKET WARRANTY (excluding lenses): 4 years against breakage if properly cared for. **Proper care entails hydrating (soaking) the nylon frame in water for 24 hours every 2-3 months.**

